

F&E MODULE

Using the ACMS™ F&E Module enables the scanning or loading of any number or type of samples. There are two evaluation options available for the selected samples: standardized evaluation and the evaluation colour comparison.

The standardized evaluation allows the user to calculate pixel-by-pixel colour values on an arbitrary ROI (region of interest). The data is saved in an exported evaluation file.

The evaluation colour comparison enables the user to compare two arbitrary ROIs from a single sample (e.g. in gradient levels). Thus, it is possible to calculate up to three comparisons simultaneously. The absolute and relative mean values are displayed for each evaluated pair. The data can be exported.

The generated ACMS™ measurement data provides the user with an objective database and hence with the possibility of carrying out further evaluations at the user's own discretion.